

Anti-MDCK HCP Antibody

Product No.: 1301308-Ab01

■Product Details

Size	1 mg
Species Reactivity	<i>Canis familiaris</i> , dog
Host/Isotype	Sheep / IgG
Class	Polyclonal
Type	Antibody
Conjugate	Unconjugated
Immunogen	MDCK HCPs produced by sonication cell lysis
Form	Liquid
Purification	Affinity Purified
Buffer	PBS
Contains	No preservatives
Applications	HCP Coverage evaluations, Western Blot, ELISA
Storage	Maintain at 2-8°C for short-term use. For long-term storage, please store below -65°C.
Validity period	12 months

■Attentions

This product is for research use only but not intended for diagnostic purposes.

■Related Product

Product Name	Product No.	Quantity
MDCK HCP ELISA Kit (One-step ELISA)	1301308	96 Tests

Effective date: 10 Jul. 2024

Support & Contact

SHENTEK

Huzhou Shenke Biotechnology Co., Ltd.
www.shentekbio.com
Address: 8th Floor, 6B Building, No.1366 Hongfeng Road, Huzhou313000, Zhejiang Province, China
E-mail: info@shentekbio.com
Phone: +1 (908) 822-3199 / (+86) 400-878-2189