

# REO-3 Positive Control

**Product No.: 1506737**

**(For Research Use Only)**

## ■ Product description

SHENTEK® REO-3 Positive Control is intended for use in combination with SHENTEK® Virus DNA & RNA Extraction and Detection Kits, and serves as a positive control to monitor the entire process of REO-3 RNA extraction and RT-PCR detection assay.

The product described herein is for research use only and is not intended for diagnostic or therapeutic use.

## ■ Kit contents and storage

Reagent	Quantity	Storage
REO-3 Positive Control	100 µL × 3 tubes	below -65°C

REO-3 Positive Control can be stored at appropriate conditions for up to 12 months.

Please check the expiration date on the labels.

## ■ Attentions

1. Thaw one tube completely at 2-8°C or on ice, then vortex thoroughly and centrifuge briefly.
2. Follow the Positive Control Sample (PCS) preparation protocol outlined in the user guide of SHENTEK® REO-3 Detection Kit or SHENTEK® Bovine Virus Detection Kit (BVDV/REO-3/PI-3/BPV-3/BAV-3).
3. This product is ready for direct use without dilution.
4. Use only after it is equilibrated to room temperature, and avoid repeated freeze-thaw cycles.
5. Do not mix reagents from different batches.
6. This product is an armored RNA, which is safe, but it is still recommended to operate it in a biosafety cabinet. Experimental waste generated during its use should be sterilized via autoclaving and disposed of in compliance with medical waste regulations.

**■ Related products**

Product Name	Product No.	Quantity
SHENTEK <sup>®</sup> Virus DNA & RNA Extraction Kit	1506730	50 Extractions
SHENTEK <sup>®</sup> REO-3 Detection Kit	1506733	50 Reactions
SHENTEK <sup>®</sup> Bovine Virus Detection Kit (BVDV/REO-3/PI-3/BPV-3/BAV-3)	1506741	50 Reactions

Effective date: 26 Feb. 2026

**Support & Contact****SHENTEK**

Huzhou Shenke Biotechnology Co., Ltd.

[www.shentekbio.com](http://www.shentekbio.com)

Address: 8th Floor, 6B Building, No.1366 Hongfeng Road, Huzhou 313000, Zhejiang Province, China

E-mail: [info@shentekbio.com](mailto:info@shentekbio.com)

Phone: +1 (908) 822-3199 / (+86) 400-878-2189