

Anti-CHO-K1 HCP Antibody

Product No.: 1301305-Ab02

■ Product information

Size	1 mg
Species Reactivity	Chinese hamster
Host/Isotype	Sheep / IgG
Class	Polyclonal
Type	Antibody
Conjugate	Unconjugated
Immunogen	CHO HCPs produced by CHO K1 fed-batch supernatant
Form	Liquid
Purification	Affinity purified
Buffer	PBS
Contains	No preservatives
Applications	HCP coverage evaluations, Western Blot, ELISA
Storage	Maintain at 2-8°C for short-term use. For long-term storage, please store below -65°C.
Validity period	12 months

■ Attentions

This product is for research use only but not intended for diagnostic purposes.

■ Related product

Product Name	Product No.	Quantity
CHO-K1 HCP ELISA Kit (One-step ELISA)	1301305-1	96 tests/Kit

Effective date: 01 Jul. 2025

Support & Contact

SHENTEK

Huzhou Shenke Biotechnology Co., Ltd.

www.shentekbio.com

Address: 8th Floor, 6B Building, No.1366 Hongfeng Road, Huzhou 313000, Zhejiang Province, China

E-mail: info@shentekbio.com

Phone: +1 (908) 822-3199 / (+86) 400-878-2189