

# Anti-*P.pastoris* HCP Antibody

Product No.: 1301313-Ab01

## ■ Product information

<b>Size</b>	1 mg
<b>Species Reactivity</b>	<i>Pichia pastoris</i> (i.e., <i>Komagataella phaffii</i> )
<b>Host/Isotype</b>	Sheep / IgG
<b>Class</b>	Polyclonal
<b>Type</b>	Antibody
<b>Conjugate</b>	Unconjugated
<b>Immunogen</b>	<i>P.pastoris</i> HCPs produced by high-pressure homogenization
<b>Form</b>	Liquid
<b>Purification</b>	Affinity purified
<b>Buffer</b>	PBS
<b>Contains</b>	No preservatives
<b>Applications</b>	HCP coverage evaluations, Western Blot and ELISA
<b>Storage</b>	Maintain at 2-8°C for short-term use. For long-term storage, please store below -65°C.
<b>Validity period</b>	12 months

## ■ Attentions

This product is for research use only but not intended for diagnostic purposes.

## ■ Related product

Product Name	Product No.	Quantity
<i>P.pastoris</i> HCP ELISA Kit (One-step ELISA)	1301313	96 Tests

Effective date: 08 Jul. 2024

## Support & Contact

**SHENTEK**

Huzhou Shenke Biotechnology Co., Ltd.

[www.shentekbio.com](http://www.shentekbio.com)

Address: 8th Floor, 6B Building, No.1366 Hongfeng Road, Huzhou313000, Zhejiang Province, China

E-mail: [info@shentekbio.com](mailto:info@shentekbio.com)

Phone: +1 (908) 822-3199 / (+86) 400-878-2189